


# AUTOHALER

An autohaler is breath-activated aerosol that delivers a measured dose of medication directly to the airways

## Recommended age of use

7 years and older

<p>1. Remove the mouthpiece cover from the autohaler.</p>	 <p>The diagram shows an autohaler device. A grey lever is located at the top. A mouthpiece cover is shown detached to the left. The autohaler has air holes at the base. The device is shown in an upright position.</p>
<p>2. Push the grey lever on top of the autohaler into an upright position and shake the autohaler.</p>	
<p>3. Breathe out slowly, away from the mouthpiece.</p>	<p>4. Place mouthpiece between lips and breathe in slowly and deeply.</p>
<p>5. Remove the autohaler from the mouth and hold your breath for about 10 seconds. Breathe out slowly through the nose.</p>	<p>6. Push the lever back down.</p> <p>*If more than one dose is required, repeat steps 2-6.</p> <p>Replace mouthpiece cover.</p>

- A clicking sound will be heard when taking a deep breath. This is the medication being delivered. It is important that you do not stop breathing in when this is heard
- Ensure that you do not cover the air holes at the base of the autohaler