

Cleaning Instructions For Spacer Devices And Puffers

Spacers

- Dismantle the spacer.
- Wash all the parts in clean, warm soapy water.
- Do not rinse.
- Allow the pieces to air dry. (They **must not** be dried with a cloth as this causes static electricity to build up, resulting in the medicine clinging to the inside wall of the spacer).
- Reassemble the pieces when dry.

It is recommended that you clean your spacer every 1-2 weeks

Puffers

- Remove the metal canister from the plastic holder.
- Do not wet metal canister or immerse in water.
- Rinse the plastic holder under warm running water.
- Completely dry plastic holder inside and out.
- Reassemble the canister in the holder.
- Keep the cap in place when not using puffer.

It is recommended that you clean your puffer weekly* and whenever it does not appear to be spraying well

Expiry date can be found on the side of the metal canister and this should be checked regularly

****Intal Forte and Tilade puffers should be cleaned daily and the manufacturers instructions followed.***